PORTUGUÊS

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
English

Aspiration Procedure:

1. Before use, check that the packaging is intact and undamaged.
2. Check that the length and diameter of the needle correspond to those indicated on the label.
3. Open the packaging and remove the protective cap.
4. Disinfect the site and administer periosteal local anesthesia at the puncture site. If necessary, use a sharp scalpel to make an incision at the puncture site to facilitate penetration through the skin.
5. Hold the needle with the proximal end in the palm and the index finger against the cannula near the tip. This position stabilizes the needle and allows for better control.
6. Insert the needle through the incision pointing towards the anterior superior iliac spine caudally at an angle of 15° or into the anterior iliac crest cranially at an angle of 15°.
7. Pressing gently but firmly, advance the tip of the needle through the periosteum, into the cortex, rotating the needle in a clockwise - counterclockwise motion.
8. Remove the stylet by turning the upper section of the handle 90° and pulling it out.
9. Continue to advance the needle cannula forward, during alternating clockwise and counterclockwise rotation, into the bone marrow cavity. Entry into the marrow cavity is usually detected by a decrease in resistance.
10. Connect a syringe with a luer lock to the needle.
11. Apply suction by pulling the plunger of the syringe. Remove the aspirated sample from the syringe.
12. If additional sample aspiration is required, insert the second stylet into the cannula and lock it with a counterclockwise rotation. Continue to insert the needle into the bone marrow cavity with an alternating movement. Remove the stylet and proceed as in step 10.
13. Remove the needle from the bone with light movements and treat the wound site.
14. Dispose of the needle in an appropriate container in accordance with specific hospital procedures.
15. If the needle is equipped with NeedleVISE®, insert the tip into the NeedleVISE® and dispose of in an appropriate container in accordance with specific hospital procedures.

Warnings and Precautions:

To be used only by specialized medical practitioners. The using procedure has an informative purpose, specialist physician is responsible for the proper procedure and techniques to be used with this device; eventual variations are at specialist physician's discretion. Verify the expire date before use. Used medical devices have a high contamination risk, be very careful to avoid accidental puncture. Treat the medical devices as contaminated hospital waste. Dispose in accordance with applicable laws and regulations.

ATTENTION THE NEEDLE IS NOT SUITABLE FOR USE UNDER MRI (Magnetic Resonance imaging)

This product has been explicitly designed and manufactured for single use only. The reuse can cause the loss of mechanical properties and invalidate device's functionality. The company M.D.L. SRL will not be held responsible in any way for the consequences in case of reuse of the device.

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